



## Single point

### 1. Set points in reference image

- select 'Points'
- scroll through slices to find matching landmark
- ensure left image is marked (white border)
- press 'Set' or S on keyboard
- set points: ENTER finalizes set points, DEL deletes set points

### 2. Set points in moving images

- scroll through slices to find matching landmark
- ensure right image is marked (white border)
- press 'Set' or S on keyboard
- set points: ENTER finalizes set points, DEL deletes set points
- press '+' button or SPACE to change moving image

### 3. Finalize points

- if same amount of point is set in all images:
- press 'Fix' or F on keyboard
  - points are fixed in table

### 4. Correct or delete points

- correct single point: mark corresponding cell in table and move point in image
- delete all points with number X: enter X into field below table and press

### 5. Save landmarks (points, lines, ROIs)

- press to save all landmarks

### 6. Compare landmarks (points, lines, ROIs)

- press 'Compare' button → new GUI will pop up with the results

#### Points

- Left-click sets point
- Right-click sets last point
- ENTER finalizes input
- DEL or BACKSPACE deletes last point
- ESC aborts action

## Line/ROI

### 1a. Set line/ROI in reference image

- select 'Lines' or 'ROIs'
- scroll through slices to find matching landmark
- ensure left image is marked (white border)
- press 'Set' or S on keyboard
- set lines/ROIs: 2x left mouse-click finalizes set line/ROI, DEL deletes complete line/ROI

### 1b. Correct or fix line/ROI in reference image

- as long as you see a yellow line, you can move the individual points
- press 'Fix' or F on keyboard to fix line/ROI

### 2a. Set line/ROI in moving images

- scroll through slices to find matching landmark
- ensure right image is marked (white border)
- press 'Set' or S on keyboard
- set lines/ROIs: 2x left mouse-click finalizes set line/ROI, DEL deletes complete line/ROI

### 2b. Correct or fix line/ROI in moving images

- as long as you see a yellow line, you can move the individual points
- press 'Fix' or F on keyboard to fix line/ROI
- press '+' button or SPACE to change moving image

### 3. Delete line/ROI in moving images

- delete all lines/ROIs with number X: enter X into field below table and press

#### Lines/ROIs

- Left-click sets new point
- Double-click sets last point and fixes input
- Right-click finalizes input (without fixing)
- A + Left-click sets adds a new point
- DEL, BACKSPACE or ESC aborts action